

## Datasheet

### DHRS2 polyclonal antibody

**Catalog Number:** PAB17690

**Regulation Status:** For research use only (RUO)

**Product Description:** Rabbit polyclonal antibody raised against synthetic peptide of DHRS2.

**Immunogen:** A synthetic peptide corresponding to internal of human DHRS2.

**Host:** Rabbit

**Reactivity:** Human, Mouse

**Applications:** ELISA, WB-Ce

(See our web site product page for detailed applications information)

**Protocols:** See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Specificity:** This antibody detects endogenous levels of total DHRS2 protein.

**Form:** Liquid

**Recommend Usage:** Western Blot (1:500-1:1000)

ELISA (1:40000)

The optimal working dilution should be determined by the end user.

**Storage Buffer:** In PBS, pH 7.4 (150mM NaCl, 0.02% sodium azide, 50% glycerol)

**Storage Instruction:** Store at -20°C.

Aliquot to avoid repeated freezing and thawing.

**Entrez GeneID:** 10202

**Gene Symbol:** DHRS2

**Gene Alias:** HEP27, SDR25C1