

Datasheet

GHITM polyclonal antibody

Catalog Number: PAB17736

Regulation Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against synthetic peptide of GHITM.

Immunogen: A synthetic peptide corresponding to internal of human GHITM.

Host: Rabbit

Reactivity: Human

Applications: ELISA, WB-Ce

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: This antibody detects endogenous levels of total GHITM protein.

Form: Liquid

Recommend Usage: Western Blot (1:500-1:1000)

ELISA (1:10000)

The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS, pH 7.4 (150mM NaCl, 0.02% sodium azide, 50% glycerol)

Storage Instruction: Store at -20°C.

Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 27069

Gene Symbol: GHITM

Gene Alias: DERP2, DKFZp566C0746, FLJ26584, HSPC282, MICS1, My021, PTD010, TMBIM5