

Datasheet

FLI1 monoclonal antibody (M07), clone 1H4

Entrez GeneID: 2313**Gene Symbol:** FLI1**Gene Alias:** EWSR2, SIC-1**Catalog Number:** H00002313-M07**Regulation Status:** For research use only (RUO)**Product Description:** Mouse monoclonal antibody
raised against a full length recombinant FLI1.**Clone Name:** 1H4**Immunogen:** FLI1 (AAH01670, 1 a.a. ~ 452 a.a)
full-length recombinant protein with GST tag. MW of the
GST tag alone is 26 KDa.**Sequence:**

MDGTIKEALSVVSDGQSLFDSAYGAAAHLPKADMAAS
GSPDYGQPHKINLPPQQEWISQPVRVNVKREYDHM
NGSRESPVDCSVSKCSKLVGGGESNPMNYSYMDEK
NGPPPPNMTTNERRVIVPADPTLWTQEHVRQWLEWA
IKEYSLMEIDTSFFQNMDGKELCKMNKEDFLRATTLYN
TEVLLSHLSYLRESSLLAYNTTSHTDQSSRLSVKEDPS
YDSVRRGAWGNMNSGLNKSPPLGGAQTISKNTEQR
PQPDYQILGPTSSRLANPGSGQIQLWQFLLELLSDSA
NASCITWEGTNGEFKMTDPDEVARRWGERISKPNMN
YDKLSRALRYYYDKNIMTKVHGKRYAYKFDHFHIAQAL
QPHPTESMYKYPDISYMPHYHAHQKQVNFVPPHPS
SMPVTSSSFFGAASQYWTSPGGIYPNPVPRHPNTH
VPSHLGSYY

Host: Mouse**Reactivity:** Human**Applications:** ELISA, IF, WB-Ce, WB-Re(See our web site product page for detailed applications
information)**Protocols:** See our web site at<http://www.abnova.com/support/protocols.asp> or product
page for detailed protocols**Isotype:** IgG2a Kappa**Storage Buffer:** In 1x PBS, pH 7.4**Storage Instruction:** Store at -20°C or lower. Aliquot to
avoid repeated freezing and thawing.