

Datasheet

FPR2 purified MaxPab mouse polyclonal antibody (B02P)

Catalog Number: H00002358-B02P

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length human FPR2 protein.

Immunogen: FPR2 (NP_001005738.1, 1 a.a. ~ 351 a.a) full-length human protein.

Sequence:

METNFSSTPLNEYEEVSYESAGYTVLRILPLVVLGVTFVL
GVLGNGLVIWVAGFRMTRTVTTICYLNLALADFSFTAT
LPFLIVSMAMGEKWPFGWFLCKLIHIVVDINLFGSVFLI
GFIALDRCICVLHPVWAQNHRTVSLAMKVIVGPWILAL
VLTLPVFLFLTTVTIPNGDTYCTFNFAWGGTPEERLK
VAITMLTARGIIRFVIGFSLPMSIVAICYGLIAAKIHKKGMI
KSSRPLRVLTAVVASFFICWFPFQLVALLGTVWLKEML
FYGKYIIDILVNPTSSLAFFNSCLNPMLYVFGQDFRE
RLIHSPLTSLERALSEDSAPTNDTAANSASPPAETELQ
AM

Host: Mouse

Reactivity: Human

Applications: Det Ab, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 2358

Gene Symbol: FPR2

Gene Alias: ALXR, FMLP-R-II, FMLPX, FPR2A, FPRH1, FPRH2, FPRL1, HM63, LXA4R