

Datasheet

GATA3 purified MaxPab rabbit polyclonal antibody (D01P)

Catalog Number: H00002625-D01P

Regulation Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against a full-length human GATA3 protein.

Immunogen: GATA3 (NP_001002295.1, 1 a.a. ~ 444 a.a) full-length human protein.

Sequence:

MEVTADQPRWVSHHHPAVLNGQHPDTHHPGLSHSY
MDAAQYPLPEEVDVLFNIDGQGNHVPPYYGNSVRATV
QRYPPTHHGSQVCRPPLLHGSLPWLDGGKALGSHHT
ASPWNLSPFSTSIHHGSPGPLSVYPPASSSSLSGGH
ASPHLFTFPPTPPKDVSPDPSLSTPGSAGSARQDEKE
CLKYQVPLPDSMKLESSHSRGSMTALGGASSSTHHPI
TTYPPYVPEYSSGLFPPSSLLGGSPGTGFGCKSRPKAR
SSTEGRECVNCGATSTPLWRRDGTGHYLCNACGLYH
KMNGQNRPLIKPKRRLSAARRAGTSCANCQTTTTTLW
RRNANGDPVCNACGLYYKLHNINRPLTMKKEGIQTRN
RKMSKSKKCKKVHDSLEDFPKNSSFNPAALSRHMSS
LSHISPFSSHMLTTPTPMHPPSSLSFGPHHPSSMV
TAMG

Host: Rabbit

Reactivity: Human

Applications: WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 2625

Gene Symbol: GATA3

Gene Alias: HDR, MGC2346, MGC5199, MGC5445