

Datasheet

CTNNAL1 MaxPab rabbit polyclonal antibody (D01)

Catalog Number: H00008727-D01

Regulation Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against a full-length human CTNNAL1 protein.

Immunogen: CTNNAL1 (NP_003789.1, 1 a.a. ~ 734 a.a) full-length human protein.

Sequence:

MAASPGPAGVGGAGAVYGS GSSGFALDSGLEIKTRS
VEQTLLPLVSQITTLINHKDNTKKSDKTLQAIQRVGQAV
NLAVGRFVKVGEAIAENENWDLKEEINIACIEAKQAGETI
AALTDITNLNHLSDGQITIFTDKTGVIAARLLLSSVTK
VLLLADRVVIKQIITSRNKVLATMERLEKVN SFQEFVQIF
SQFGNEMVEFAHLSGDRQNDLKDEKKKAKMAAARAV
LEKCTMMLLTASKTCLRHPNCE SAHKNKEGVFDRMKV
ALDKVIEIVTDCKPNGETDISSISIFTGIKEFKMNIEALRE
NLYFQSKENLSVTLEVILERMEDFTDSAYTSHEHRERIL
ELSTQARMELQQLISVWIAQSKKTKSIAEELELSILKIS
HSLNELKKELHSTATQLAADLLKYHADHVVLKALKLTG
VEGNLEALAEYACKLSEQKEQLVETCRLLRHISGTEPL
EITCIHAEETFQVTGQQIISAAETLTLHPSSKIAKENLDV
FCEAWESQISDMSTLLREINDVFEGRRGEKYGYLSLP
KPMKNNANLKS LKPDKPDSEEQAKIAKLGLKLGLTSD
ADCEIEKWEDQENEIVQYGRNMSSMAYSLYLFTRGEG
PLKTSQDLIHQLEVFAAEG LKTSSVQAFSKQLKDDDK
LMLLLEINKLIPLCHQLQTVTKTSLQNKVFLKVDKCITKT
RSMALLVQLLSLCYKLLKKLQ MENNGWVSVTNKDT
MDSKT

Host: Rabbit

Reactivity: Human

Applications: IP, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: No additive

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 8727

Gene Symbol: CTNNAL1

Gene Alias: CLLP, FLJ08121, alpha-CATU