

Datasheet

N4BP1 purified MaxPab mouse polyclonal antibody (B01P)

Catalog Number: H00009683-B01P

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length human N4BP1 protein.

Immunogen: N4BP1 (AAI60019.1, 1 a.a. ~ 896 a.a) full-length human protein.

Sequence:

MAARAVLDEFTAPAEKAEELLEQSRGRIEGLFGVSLAVL
GALGAEELPARIWLQLCGAQEAHVSAKEYIKGICEPE
LEERECYPKDMHCIFVGAESLFLKSLIQDTCADLCILDI
GLLGIRGSAEAVVMARSHIQFVKLFENKENLPSSQKE
SEVKREFKQFVEAHADNYTMDLLIPTSLKKELLTLTQ
GEENLFETGDDEVIMRDSQQTEFTQNAATGLNISRD
ETVLQEEARNKAGTPVSELTQMDTVLSSSPDVLFDPI
NGLTPDEEALSNERICQKRRFSDSEERHTKKQFSLEN
VQEGEILHDAKTLAGNVIADLSDSSADSENLSPIKET
EEMEYNILVNFFKTMGYEQIVEKVIKVGYPSTEPLLLL
EEIEKENKRFQEDREFSAGTVYPETNKTGNKGVYSST
NELTTDSTPKKTQAHTQQNMVEKFSQLPFKVEAKPCT
SNCRINTFRTPVIEQKHEVWGSNQNYICNTDPETDGLS
PSVASPSPKEVNFVSRGASSHQPRVPLFPENGLHQQ
PEPLLNNMKSACEKRLGCCSSPHSKPNCSTLSPPMP
LPQLLPSVTDARSAGPSDHIDSSVTGVQRFRTDLKIPY
KLELKNEPGRDLDKHIVIDGSNVAITHGLKKFFSCRGIAI
AVEYFWKLGNRNITVFVPQWRTRRDPNVTEQHFLTQL
QELGILSLTPARMVFGERIASHDDRFLHLADKTGGIIV
TNDNFREFVNESVSWREIITKRLQYTFVGDIFMVPDD
PLGRSGPRLEEFLLQKEVCLRDMQPLLSALPNVGMFDP
SFRVPGTQAASSTSHQPPTRIQGAPSSHWLPQQPHFPL
LPALPSLQQNLPMQAQRSSAETNELREALLKIFPDSEQ
RLKIDQILVAHPYMKDLNALSAMVLD

Host: Mouse

Reactivity: Human

Applications: Det Ab, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product

page for detailed protocols

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 9683

Gene Symbol: N4BP1

Gene Alias: FLJ31821, KIAA0615, MGC176730