

Datasheet

SPATA2 MaxPab rabbit polyclonal antibody (D01)

Gene Alias: FLJ13167, KIAA0757, PD1, tamo

Catalog Number: H00009825-D01

Regulation Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against a full-length human SPATA2 protein.

Immunogen: SPATA2 (NP_006029.1, 1 a.a. ~ 520 a.a) full-length human protein.

Sequence:

MGKPSSMDTKFKDDLFRKYVQFHESKVDTTTSRQRP
GSDECLRVAASTLLSLHKVDPFYRFRLIQFYEVVSSL
RSLSSSLRALHGAFSMLTVGINFLYPWKKEFRSIK
TYTGPVYVVKSTLLEEDIRAILSCMGYTPELGTAYKLR
ELVETLQVKMVSFELFLAKVECEQMLEIHSQVKDKGYS
ELDIVSERKSSAEDVRGCSDALRRRAEGREHLTASMS
RVALQKSASERAADYKPRVTKPSRSVDAYDSYWE
SRKPPLKASLSLRKEPVATDVGDDLKDEIRPSPSLLTM
ASSPHGSPDVLPPASPSNGPALLRGTYFSTQDDVDLY
TDSEPRATYRRQDALRPDVLRLNDAHSLYHKRSPP
AKESALSCKQSCGLSCSSSLCQRCDSLLTCPPASKPS
AFPSKASTHDSLHAGASLREKYPGQTQGLDRLPHLHS
KSKPSTTPTSRCGFCNRPGATNTCTQCSKVSCDACLS
AYHYDPCYKKSELHKFMPNNQLNYKSTQLSHLVYR

Host: Rabbit

Reactivity: Human

Applications: IP, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: No additive

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 9825

Gene Symbol: SPATA2