

Datasheet

PSMD6 purified MaxPab mouse polyclonal antibody (B01P)

Catalog Number: H00009861-B01P

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length human PSMD6 protein.

Immunogen: PSMD6 (NP_055629.1, 1 a.a. ~ 389 a.a) full-length human protein.

Sequence:

MPLENLEEEGLPKNPDLRIAQLRFLSLPEHRGDAAVR
DELMAAVRDNNMAPYYEALCKSLDWQIDVDLLNKMKK
ANEDELKRLDEELEDAAEKNLGESEIRDAMMAKAEYLC
RIGDKEGALTAFRKTYDKTVALGHRDLIVFYLLRIGLFY
MDNDLITRNTEKAKSLIEEGGDWDRRNRLKVYQGLYC
VAIRDFKQAAELFLDTVSTFTSYELMDYKTFVTYTVYV
SMIALERPDLREKVIKGAIEVLHSLPAVRQYLFSLYE
CRYSVFFQSLAVVEQEMKKDWLFAPHYRYVREMRIH
AYSQLESYRSLTLGYMAEAFGVGVFEFIDQELSRFIAA
GRLHCKIDKVNEIVETNRPDSKNWQYQETIKKGDLLN
RVQKLSRVINM

Host: Mouse

Reactivity: Human

Applications: WB-Ti, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 9861

Gene Symbol: PSMD6

Gene Alias: KIAA0107, Rpn7, S10, SGA-113M, p44S10