

Datasheet

ANGPTL7 monoclonal antibody (M02), clone 3E7

Catalog Number: H00010218-M02

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody
raised against a partial recombinant ANGPTL7.

Clone Name: 3E7

Immunogen: ANGPTL7 (NP_066969, 27 a.a. ~ 126 a.a)
partial recombinant protein with GST tag. MW of the
GST tag alone is 26 KDa.

Sequence:

QKLSKHKTPAQQLKAANCCEEVKELKAQVANLSSLL
SELNKKQERDWVSVVMQVMELESNSKRME SRLTDAE
SKYSEMNNQIDIMQLQAAQTVTQTSAD

Host: Mouse

Reactivity: Human

Applications: ELISA, IHC-P, WB-Re
(See our web site product page for detailed applications
information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product
page for detailed protocols

Isotype: IgG2a Kappa

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to
avoid repeated freezing and thawing.

Entrez GeneID: 10218

Gene Symbol: ANGPTL7

Gene Alias: AngX, CDT6, RP4-647M16.2, dJ647M16.1