

Datasheet

ANP32B monoclonal antibody (M05), clone 8E10

Catalog Number: H00010541-M05

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody
raised against a full-length recombinant ANP32B.

Clone Name: 8E10

Immunogen: ANP32B (NP_006392.1, 1 a.a. ~ 251 a.a)
full-length recombinant protein with GST tag. MW of the
GST tag alone is 26 KDa.

Sequence:

MDMKRRIHLELRNRTPAAVRELVLDNCKSNDGKIEGLT
AEFVNLEFLSLINVGLISVSNLPKLPKPKKLELSENRIFG
GLDMLAEKLPNLTHLNLSGNKLDISTLEPLKKLECLKS
LDLFNCEVTNLNDYRESVFKLLPQLTYLDGYDREDQE
APDSDAEVDGVDEEEEEDEGEDEDEDEDEDEGEDEEEF
DEEDEDEDEVEGDEDDDEVSEEEEEFGLDEEDEDED
EDEEEEEGGKGEKRETDDEGEDD

Host: Mouse

Reactivity: Human

Applications: ELISA, IF, S-ELISA, WB-Re, WB-Tr
(See our web site product page for detailed applications
information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product
page for detailed protocols

Isotype: IgG2a Kappa

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to
avoid repeated freezing and thawing.

Entrez GeneID: 10541

Gene Symbol: ANP32B

Gene Alias: APRIL, PHAPI2, SSP29