

Datasheet

IVNS1ABP monoclonal antibody (M01), clone 4B1

Catalog Number: H00010625-M01

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against a partial recombinant IVNS1ABP.

Clone Name: 4B1

Immunogen: IVNS1ABP (NP_006460.2, 1 a.a. ~ 100 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.

Sequence:

MIPNGYLMFEDENFIESSVAKLNALRKSGQFCDVRLQV
CGHEMLAHRAVLACCSPYLFEIFNSDSDPHGISHVKFD
DLNPEAVEVLLNYAYTAQLKADKE

Host: Mouse

Reactivity: Human

Applications: ELISA, S-ELISA, WB-Re
(See our web site product page for detailed applications information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Isotype: IgG1 Kappa

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 10625

Gene Symbol: IVNS1ABP

Gene Alias: DKFZp686K06216, FLARA3, FLJ10069, FLJ10411, FLJ10962, FLJ35593, FLJ36593, HSPC068, KIAA0850, ND1, NS-1, NS1-BP, NS1BP