

Datasheet

MALT1 MaxPab rabbit polyclonal antibody (D01)

Catalog Number: H00010892-D01

Regulation Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against a full-length human MALT1 protein.

Immunogen: MALT1 (NP_776216.1, 1 a.a. ~ 813 a.a) full-length human protein.

Sequence:

MSLLGDPLQALPPSAAPTGPLLAPPAGATLNRLREPLL
RRLSELDDQAEGRGWRRRLAELAGSRGRLRLSCLDLE
QCSLKVLEPEGSPSLCLLKLMEKGCTVTELSDFLQA
MEHTEVLQLLSPPGIKITVNPESKAVLAGQFVKLCRA
TGHPFVQYQWFKMNKEIPNGNTSELIFNAVHVKDAGF
YVCRVNNNFTFEFSQWSQLDVCDIPESFQRSVDGVSE
SKLQICVEPTSQKLMPGSTLVLCVAVGSPIPHYQWF
KNELPLTHETKKLYMVPYVDLEHQGTYWCHVYNDRD
SQDSKKVEIIDEELNNLGHDPNKEQTTDQPLAKDKVALL
IGNMNYREHPKLKAPLVDVYELTNLLRQLDFKVVSLLD
LTEYEMRNAVDEFLLLDKGVYGLLYAGHGYENFGN
SFMVPVDAPNPYRSENCLCVQNILKLMQEKETGLNVF
LLDMCRKRNDYDDTIPILDALKVTANIVFGYATCQGA
AFEIQHSLGLANGIFMKFLKDRILLEDKKITVLLDEVAEDM
GKCHLTGKGQALEIRSSLSEKRALTDPIQGTEYSAESL
VRNLQWAKAHELPEMCKLFCGVCQIQLGFAAEFSNV
MIIYTSIVYKPEIIMCDAYVTFPLDLIDPKDANKGTP
EETGSYLVS KDLPKHCLYTRLSSLQKLKEHLVFTVCLS
YQYSGLEDTVEDKQEVNVGKPLIAKLDMHRGLGRKTC
FQTCLMSNGPYQSSAATSGGAGHYHSLQDPFHGVYH
SHPGNPSNVTPADSCHCSRTPDAFISSFAHHASCHFS
RSNVPVETTDEIPFSFSDRLRISEK

Host: Rabbit

Reactivity: Human

Applications: IF, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: No additive

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 10892

Gene Symbol: MALT1

Gene Alias: DKFZp434L132, MLT, MLT1

Gene Summary: This gene has been found to be recurrently rearranged in chromosomal translocation with two other genes - baculoviral IAP repeat-containing protein 3 (also known as apoptosis inhibitor 2) and immunoglobulin heavy chain locus - in mucosa-associated lymphoid tissue lymphomas. The protein encoded by this gene may play a role in NF-kappaB activation. Two alternatively spliced transcript variants encoding different isoforms have been described for this gene. [provided by RefSeq]