

Datasheet

MFGE8 polyclonal antibody (A01)

Catalog Number: H00004240-A01

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a partial recombinant MFGE8.

Immunogen: MFGE8 (NP_005919, 61 a.a. ~ 170 a.a) partial recombinant protein with GST tag.

Sequence:

YAGNHCETKCVPLGMENGNIANSQIAASSVRVTLG
LQHWVPELARLNLRAGMVNAWTPSSNDDNPWIQVNLL
RRMWVTGVVTQGASRLASHEYLKAFKVAYSLNGHEF
D

Host: Mouse

Reactivity: Human

Applications: ELISA, WB-Re

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: 50 % glycerol

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 4240

Gene Symbol: MFGE8

Gene Alias: BA46, EDIL1, HsT19888, OAcGD3S, SED1, SPAG10, hP47