

Datasheet

MLF1 purified MaxPab rabbit polyclonal antibody (D01P)

Catalog Number: H00004291-D01P

Regulation Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against a full-length human MLF1 protein.

Immunogen: MLF1 (NP_071888.1, 1 a.a. ~ 268 a.a) full-length human protein.

Sequence:

MFRMLNSSFEDDPFFSESILAHRENMRQMIRSFSEPF
GRDLLSISDGRGRAHNRRGHNDGEDSLTHTDVSSFQ
TMDQMVSNMRNYMQKLERNFGQLSVDPNGHSFCSS
SVMTYSKIGDEPPKVFQASTQTRRAPGGIKETRKAMR
DSDSGLEKMAIGHHHIHDRAHVIKKSKNKTGDDEEVNQ
EFINMNESDAHAFDEEWQSEVLKYKPGRHNLGNTRM
RSVGHENPGSRELKRREKPQQSPAIEHGRRSNVLGD
KLHIKGSSVKSNNK

Host: Rabbit

Reactivity: Human, Mouse

Applications: WB-Ti, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 4291

Gene Symbol: MLF1

Gene Alias: -