

Datasheet

MAP3K9 polyclonal antibody (A01)

Catalog Number: H00004293-A01

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a partial recombinant MAP3K9.

Immunogen: MAP3K9 (NP_149132, 971 a.a. ~ 1080 a.a) partial recombinant protein with GST tag.

Sequence:

EFRLPDPNVVFPPTPRRWNTQQDSTLERPKTLEFLP
RPRPSANRQLDPWWFVSPSHARSTSPANSSSTETP
SNLDS CFASSSSTVEERPGLPALLPFQAGLPPTERT

Host: Mouse

Reactivity: Human

Applications: ELISA, WB-Re

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: 50 % glycerol

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 4293

Gene Symbol: MAP3K9

Gene Alias: MLK1, PRKE1