

Datasheet

MST1 MaxPab rabbit polyclonal antibody (D01)

Catalog Number: H00004485-D01

Regulation Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against a full-length human MST1 protein.

Immunogen: MST1 (NP_066278.2, 1 a.a. ~ 711 a.a) full-length human protein.

Sequence:

MGWLPLLLLLTQCLGVPGQRSPLNDFQVLRGTELQHL
LHAVVPGPWQEDVADAEACAGRCGPLMDCRAFHYN
VSSHGCQLLPWTQHSPHTRLRRSGRCDLFQKKDYVR
TCIMNNGVG YRGTMATTVGGLPCQAWSHKFPNDHKY
TPTLRNGLEENFCRNPDPGGPWCTTDPVAVRFQS
CGIKSCREAAACVWCNGEYRGAVDRTESGRECQRW
DLQHPHQHPFEPGKFLDQGLDDNYCRNPDGSERPW
CYTTDPQIEREFCDLPRCGSEAQPRQEATTVSCFRGK
GEGYRGTTANTTTAGVPCQRWDAQIPHQHRFTPEKYA
CKDLRENFCRNPDGSEAPWCFTLRPGMRAAFYQIR
RCTDDVRPQDCYHGAGEQYRGTVSKTRKGVQCQRW
SAETPHKPQFTFTSEPHAQLEENFCRNPDGDSHGPW
CYTMDPRTPFDYCALRRCADDPSPSILDPDPDQVQFEK
CGKRVDRLDQRRSKLRVVGHPGNSPWTVSLRNQ
GQHFCGSLVKEQWILTARQCFSSCHMPLTG YEVWL
GTLFQNPQHGEPSLQRPVAKMVCGPSGSQLVLLKLE
RSVTLNQRVALICLPPEWYVPPGKCEIAGWGETKG
TGNDTVLNVALLNVISNQEKNIKHRGRVRESEMCTEG
LLAPVGACEGDYGGPLACFTHNCWVLEGIIPNRVCAR
SRWPAVFTRVSVFVDWIHKVMRLG

Host: Rabbit

Reactivity: Human

Applications: IP, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: No additive

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 4485

Gene Symbol: MST1

Gene Alias: D3F15S2, DNF15S2, HGFL, MSP, NF15S2