

Datasheet

MST1R monoclonal antibody (M07), clone 1H7

Catalog Number: H00004486-M07

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody
raised against a partial recombinant MST1R.

Clone Name: 1H7

Immunogen: MST1R (NP_002438, 26 a.a. ~ 125 a.a)
partial recombinant protein with GST tag. MW of the
GST tag alone is 26 KDa.

Sequence:

DWQCPRTPTYAASRDFDVKYVVPFSAGGLVQAMVTY
EGDRNESAVFVAIRNRLHVLGPDLSVQSLATGPAGD
PGCQTCAACGPGPHGPPGDTDTKVLVL

Host: Mouse

Reactivity: Human

Applications: ELISA, S-ELISA

(See our web site product page for detailed applications
information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product
page for detailed protocols

Isotype: IgG2a Kappa

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to
avoid repeated freezing and thawing.

Entrez GeneID: 4486

Gene Symbol: MST1R

Gene Alias: CD136, CDw136, PTK8, RON