

Datasheet

MYBL1 monoclonal antibody (M01), clone 2A2

Catalog Number: H00004603-M01

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody
raised against a partial recombinant MYBL1.

Clone Name: 2A2

Immunogen: MYBL1 (XP_034274, 625 a.a. ~ 734 a.a)
partial recombinant protein with GST tag. MW of the
GST tag alone is 26 KDa.

Sequence:

KSLVLDNWEKEESGTQLLTEDISDMQSENRFSTSLMI
PLLEIHDNRCNLIPEKQDINSTNKTYTLTKKKPNPNTSK
VVKLEKNLQSNCEWETVVYGTKEDQLIMTEQAR

Host: Mouse

Reactivity: Human

Applications: ELISA, S-ELISA, WB-Ce, WB-Re
(See our web site product page for detailed applications
information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product
page for detailed protocols

Isotype: IgG1 Kappa

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to
avoid repeated freezing and thawing.

Entrez GeneID: 4603

Gene Symbol: MYBL1

Gene Alias: A-MYB, AMYB, MGC120059, MGC120061