

Datasheet

MYF6 monoclonal antibody (M05), clone 1H3

Catalog Number: H00004618-M05

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody
raised against a partial recombinant MYF6.

Clone Name: 1H3

Immunogen: MYF6 (NP_002460.1, 1 a.a. ~ 99 a.a)
partial recombinant protein with GST tag. MW of the
GST tag alone is 26 KDa.

Sequence:

MMMDLFETGSYFFYLDGENVTLQPLEVAEGSPLYPGS
DGTLSPCQDQMPPEAGSDSSGEEHVLAPPGLQPPHC
PGQCLIWACKTCKRKSAPTDRRKAAT

Host: Mouse

Reactivity: Human

Applications: ELISA, S-ELISA, WB-Re
(See our web site product page for detailed applications
information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product
page for detailed protocols

Isotype: IgG2a Kappa

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to
avoid repeated freezing and thawing.

Entrez GeneID: 4618

Gene Symbol: MYF6

Gene Alias: MRF4, bHLHc4