

Datasheet

NAP1L1 MaxPab mouse polyclonal antibody (B01)

Catalog Number: H00004673-B01

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length human NAP1L1 protein.

Immunogen: NAP1L1 (AAH02387, 1 a.a. ~ 391 a.a) full-length human protein.

Sequence:

MADIDNKEQSELDQDLDDVEEVEEETGEETKLKARQ
LTVQMMQNPQILAALQERLDGLVETPTGYIESLPRVVK
RRVNALKNLQVKCAQIEAKFYEEVHDLERKYAVLYQPL
FDKRFEIINAIYEPTEECEWKPDDEDEISEELKEKAKIE
DEKKDEEKEDPKGIPEFWLTVFKNVDLLSDMVQEHDE
PILKHLKDIKVKFSDAGQPMSFVLEFHFEPNEYFTNEVL
TKTYRMRSEPDDSDPFSFDGPEIMGCTGCQIDWKKG
KNVTLKTIKKKQKHKGRTVTVTKTVSNDFFNFAP
PEVPESGDLDDDAEAILAADFEIGHFLRERIIPRSVLYFT
GEAIEDDDDDYDEEGEEADEEGEEEGDEENDPDYDP
KKDQNP AECKQQ

Host: Mouse

Reactivity: Human

Applications: Det Ab, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: No additive

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 4673

Gene Symbol: NAP1L1

Gene Alias: FLJ16112, MGC23410, MGC8688, NAP1, NAP1L, NRP

Gene Summary: This gene encodes a member of the nucleosome assembly protein (NAP) family. This protein participates in DNA replication and may play a role in modulating chromatin formation and contribute to the regulation of cell proliferation. Alternative splicing of this gene results in several transcript variants; however, not all have been fully described. [provided by RefSeq]