

Datasheet

NAP1L2 MaxPab rabbit polyclonal antibody (D01)

Catalog Number: H00004674-D01

Regulation Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against a full-length human NAP1L2 protein.

Immunogen: NAP1L2 (AAH26325.1, 1 a.a. ~ 460 a.a) full-length human protein.

Sequence:

MAKSENKELSESSQEEAGNQIMVEGLGEHLERGED
AAAGLGDDGKCGEEAAAGLGEEGENGEDTAAGSGED
GKKGGDTDEDSEADRPKGLIGYVLDTFVESLPVKVK
YRVLALKKLQTRAANLESKFLREFHDIERKFAEMYQPL
LEKRRQIINAIYEPTEEECYKSDSEDCDDEEMCHEEM
YGNEEGMVHEYVDEDDGYEDYYDYAVEEEEEEEEE
DDIEATGEENKEEEDPKGIPDFWLTVLKNVDLTLPLIKK
YDEPILKLLTDIKVKLSDPGEPLSFTLEFHFKNPFYFKN
ELLTKTYVLKSKLAYYDPHPYRGTAIEYSTGCEIDWNE
GKNVTLTKTIKKKQKHRIWGTIRTVTEDFPKDSFFNFFSP
HGITSNGRDGNDFFLLGHNLRTYIIPRSVLFFSGDALE
SQQEGVVREVNDIAIDKIIDNWMMAAIEEVKACCKNLE
ALVEDIDR

Host: Rabbit

Reactivity: Human

Applications: IP, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: No additive

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 4674

Gene Symbol: NAP1L2

Gene Alias: BPX, MGC26243

Gene Summary: This gene encodes a member of the nucleosome assembly protein (NAP) family. The function of this family member is unknown; however, mouse studies suggest that it represents a class of tissue-specific factors interacting with chromatin to regulate neuronal cell proliferation. [provided by RefSeq]