

Datasheet

NDUFB10 polyclonal antibody (A01)

Catalog Number: H00004716-A01

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a partial recombinant NDUFB10.

Immunogen: NDUFB10 (NP_004539, 73 a.a. ~ 172 a.a) partial recombinant protein with GST tag.

Sequence:

EEDIMCMYEAEMQWKRDYKVDQEIIINIMQDRLKACQQ
REGQNYQQNCIKEVEQFTQVAKAYQDRYQDLGAYSS
ARKCLAKQRQRMLQERKAAKEAAAATS

Host: Mouse

Reactivity: Human

Applications: ELISA, WB-Ce, WB-Re

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: 50 % glycerol

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 4716

Gene Symbol: NDUFB10

Gene Alias: PDSW