

## Datasheet

### **IgG (Mouse) Isotype Control (R-PE-TXRD)**

**Catalog Number:** PAB19461

**Regulation Status:** For research use only (RUO)

**Product Description:** Mouse IgG raised from normal mouse serum.

**Host:** Mouse

**Applications:** ELISA, Flow Cyt  
(See our web site product page for detailed applications information)

**Protocols:** See our web site at  
<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Form:** Liquid

**Conjugation:** R-PE-TXRD

**Concentration:** 0.1 mg/mL

**Isotype:** IgG

**Recommend Usage:** Flow Cytometry (? 1 ug/10<sup>6</sup> cells)  
ELISA (1:500-1:5000)  
The optimal working dilution should be determined by the end user.

**Storage Buffer:** In PBS (0.09% sodium azide)

**Storage Instruction:** Store at 4°C. Do not freeze.