

Datasheet

THYN1 polyclonal antibody

Catalog Number: PAB23340

Regulation Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against recombinant THYN1.

Immunogen: Recombinant protein corresponding to amino acids of human THYN1.

Sequence:

AFFYHSNCKEPIAGLMKIVKEAYPDHTQFEKNNPHY
DPSSKEDNPKWSMVDVQFVRMMKRFIPLAELKSYHQ
AHKATGGPL

Host: Rabbit

Reactivity: Human

Applications: IF, IHC-P, WB

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Purification: Antigen affinity purification

Isotype: IgG

Recommend Usage: Immunohistochemistry

(1:200-1:500)

Western Blot (1:250-1:500)

Immunofluorescence (1-4 ug/mL)

The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS, pH 7.5 (40% glycerol, 0.02% sodium azide)

Storage Instruction: Store at 4°C. For long term storage store at -20°C.

Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 29087

Gene Symbol: THYN1

Gene Alias: HSPC144, MDS012, MGC12187, MY105, THY28, THY28KD

Gene Summary: This gene encodes a protein that is highly conserved among vertebrates and plant species and may be involved in the induction of apoptosis. Alternatively spliced transcript variants encoding different isoforms have been described. [provided by RefSeq]