



# MATERIAL SAFETY DATA SHEET

according to the Hazard Communication Standard (29 CFR 1910.1200)

Date of issue: 12/20/2012

Version 1.0

## SECTION 1. Identification

### Product identifier

Product number	818737
Product name	Chlorotrimethylsilane for synthesis

### Relevant identified uses of the substance or mixture and uses advised against

Identified uses	Chemical for synthesis
-----------------	------------------------

### Details of the supplier of the safety data sheet

Company	EMD Millipore Corporation   290 Concord Road, Billerica, MA 01821, United States of America   SDS Phone Support: +1-978-715-1335   General Inquiries: +1-978-715-4321   Monday to Friday, 9:00 AM to 4:00 PM Eastern Time (GMT-5)  e-mail: mm_sds@merckgroup.com
---------	--

Emergency telephone	800-424-9300 CHEMTREC (USA) +1-703-527-3887 CHEMTREC (International) 24 Hours/day; 7 Days/week
---------------------	--

## SECTION 2. Hazards identification

### GHS Classification

Flammable liquid, Category 2, H225  
Specific target organ systemic toxicity - single exposure, Category 3, Respiratory system, H335  
Skin corrosion, Category 1A, H314

For the full text of the H-Statements mentioned in this Section, see Section 16.

### GHS-Labeling

#### Hazard pictograms



Signal Word  
Danger

#### Hazard Statements

H225 Highly flammable liquid and vapor.  
H314 Causes severe skin burns and eye damage.

# MATERIAL SAFETY DATA SHEET

according to the Hazard Communication Standard (29 CFR 1910.1200)

Product number 818737  
Product name Chlorotrimethylsilane for synthesis

---

Version 1.0

H335 May cause respiratory irritation.

## *Precautionary Statements*

P210 Keep away from heat/sparks/open flames/hot surfaces. - No smoking.

P233 Keep container tightly closed.

P240 Ground/bond container and receiving equipment.

P280 Wear protective gloves/ protective clothing/ eye protection/ face protection.

P301 + P330 + P331 IF SWALLOWED: rinse mouth. Do NOT induce vomiting.

P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P309 + P310 IF exposed or if you feel unwell: Immediately call a POISON CENTER or doctor/physician.

P403 + P235 Store in a well-ventilated place. Keep cool.

## OSHA Hazards

This material is considered hazardous by the OSHA Hazard Communication Standard (29 CFR 1910.1200).

## Other hazards

None known.

---

## SECTION 3. Composition/information on ingredients

Formula	(CH <sub>3</sub> ) <sub>3</sub> SiCl	C <sub>3</sub> H <sub>9</sub> ClSi (Hill)
CAS-No.	75-77-4	
Molar mass	108.64 g/mol	

## Hazardous ingredients

*Chemical Name ( Concentration)*

CAS-No.

*chlorotrimethylsilane ( >= 90 % - <= 100 % )*  
75-77-4

---

## SECTION 4. First aid measures

### Description of first-aid measures

*General advice*

First aider needs to protect himself.

*Inhalation*

After inhalation: fresh air. Call in physician.

*Skin contact*

After skin contact: wash off with plenty of water. Immediately remove contaminated clothing. If available swab with polyethylene glycol 400. Call a physician immediately.

*Eye contact*

After eye contact: rinse out with plenty of water. Immediately call in ophthalmologist.

---

## MATERIAL SAFETY DATA SHEET

according to the Hazard Communication Standard (29 CFR 1910.1200)

Product number 818737  
Product name Chlorotrimethylsilane for synthesis

---

Version 1.0

### *Ingestion*

After swallowing: make victim drink water (two glasses at most), avoid vomiting (risk of perforation!). Do not attempt to neutralize. Pulmonary failure possible after aspiration of vomit. Call a physician immediately.

Never give anything by mouth to an unconscious person.

### **Most important symptoms and effects, both acute and delayed**

Irritation and corrosion, Cough, Shortness of breath, Risk of blindness!

### **Indication of any immediate medical attention and special treatment needed**

No information available.

---

## **SECTION 5. Fire-fighting measures**

### **Extinguishing media**

#### *Suitable extinguishing media*

Carbon dioxide (CO<sub>2</sub>), Dry powder

#### *Unsuitable extinguishing media*

Water, Foam

### **Special hazards arising from the substance or mixture**

Combustible.

Vapors are heavier than air and may spread along floors.

Forms explosive mixtures with air at ambient temperatures.

Pay attention to flashback.

Development of hazardous combustion gases or vapors possible in the event of fire.

Fire may cause evolution of:

Hydrogen chloride gas, silanes

Caution! in contact with water product releases:

Hydrogen chloride gas

### **Advice for firefighters**

#### *Special protective equipment for fire-fighters*

Stay in danger area only with self-contained breathing apparatus. Prevent skin contact by keeping a safe distance or by wearing suitable protective clothing.

#### *Further information*

Cool closed containers exposed to fire with water spray. Suppress (knock down) gases/vapors/mists with a water spray jet. Prevent fire extinguishing water from contaminating surface water or the ground water system.

---

## **SECTION 6. Accidental release measures**

### **Personal precautions, protective equipment and emergency procedures**

Advice for non-emergency personnel: Avoid substance contact. Do not breathe vapors, aerosols. Ensure adequate ventilation. Keep away from heat and sources of ignition. Evacuate the danger area, observe emergency procedures, consult an expert.

Advice for emergency responders: Protective equipment see section 8.

### **Environmental precautions**

Do not empty into drains. Risk of explosion.

---

# MATERIAL SAFETY DATA SHEET

according to the Hazard Communication Standard (29 CFR 1910.1200)

Product number 818737  
Product name Chlorotrimethylsilane for synthesis

---

Version 1.0

## Methods and materials for containment and cleaning up

Cover drains. Collect, bind, and pump off spills.

Observe possible material restrictions (see sections 7 and 10).

Take up with liquid-absorbent material (e.g. Chemisorb®). Dispose of properly. Clean up affected area.

---

## SECTION 7. Handling and storage

### Precautions for safe handling

Observe label precautions.

Keep workplace dry. Do not allow product to come into contact with water.

#### *Advice on protection against fire and explosion*

Keep away from open flames, hot surfaces and sources of ignition. Take precautionary measures against static discharge.

### Conditions for safe storage, including any incompatibilities

Keep away from heat and sources of ignition. Keep container tightly closed in a dry and well-ventilated place.

Store at +15°C to +25°C (+59°F to +77°F).

---

## SECTION 8. Exposure controls/personal protection

### Exposure limit(s)

Contains no substances with occupational exposure limit values.

### Engineering measures

Technical measures and appropriate working operations should be given priority over the use of personal protective equipment.

### Individual protection measures

Protective clothing should be selected specifically for the workplace, depending on concentration and quantity of the hazardous substances handled. The chemical resistance of the protective equipment should be inquired at the respective supplier.

#### *Hygiene measures*

Immediately change contaminated clothing. Apply skin- protective barrier cream. Wash hands and face after working with substance.

#### *Eye/face protection*

Tightly fitting safety goggles

#### *Hand protection*

Chemical-resistant, impervious gloves complying with an approved standard should be worn at all times when handling chemical products if a risk assessment indicates this is necessary.

#### *Other protective equipment:*

protective clothing

#### *Respiratory protection*

required when vapors/aerosols are generated.

Use a properly fitted, air-purifying or air-fed respirator complying with an approved standard if a risk assessment indicates this is necessary. Respirator selection must be based on known or anticipated exposure levels, the hazards of the product and the safe working limits of the selected respirator.

MATERIAL SAFETY DATA SHEET  
according to the Hazard Communication Standard (29 CFR 1910.1200)

Product number 818737  
Product name Chlorotrimethylsilane for synthesis

---

Version 1.0

---

**SECTION 9. Physical and chemical properties**

Physical state	liquid
Color	colorless
Odor	characteristic
Odor Threshold	No information available.
pH	No information available.
Melting point	-72 °F ( -58 °C)
Boiling point/boiling range	135 °F ( 57 °C) at 1,013 hPa
Flash point	-4 °F ( -20 °C) Method: c.c.
Evaporation rate	No information available.
Flammability (solid, gas)	No information available.
Lower explosion limit	1.5 %(V)
Upper explosion limit	46 %(V)
Vapor pressure	250 hPa at 68 °F ( 20 °C)
Relative vapor density	3.8
Relative density	0.86 g/cm <sup>3</sup> at 68 °F ( 20 °C)
Water solubility	at 68 °F ( 20 °C) (rigorous decomposition)
Partition coefficient: n-octanol/water	No information available.
Autoignition temperature	No information available.
Decomposition temperature	No information available.
Viscosity, dynamic	0.4 mPa.s at 68 °F ( 20 °C)
Explosive properties	Not classified as explosive.

## MATERIAL SAFETY DATA SHEET

according to the Hazard Communication Standard (29 CFR 1910.1200)

Product number 818737  
Product name Chlorotrimethylsilane for synthesis

---

Version 1.0

Ignition temperature 788 °F ( 420 °C)  
Method: DIN 51794

---

### SECTION 10. Stability and reactivity

#### Reactivity

Vapors may form explosive mixture with air.

#### Chemical stability

Sensitive to air.

#### Possibility of hazardous reactions

Violent reactions possible with:

Generates dangerous gases or fumes in contact with:

Water, Alcohols, Acetone, Esters, Ketones, Aldehydes, Strong acids, Bases, Amines

#### Conditions to avoid

Warming.

A range from approx. 15 Kelvin below the flash point is to be rated as critical.

Exposure to moisture.

#### Incompatible materials

no information available

#### Hazardous decomposition products

in the event of fire: See section 5.

---

### SECTION 11. Toxicological information

#### Information on toxicological effects

##### *Likely route of exposure*

Inhalation, Eye contact, Skin contact

##### *Acute oral toxicity*

LD50 rat: 4,868 mg/kg (RTECS)

Symptoms: If ingested, severe burns of the mouth and throat, as well as a danger of perforation of the esophagus and the stomach., Risk of aspiration upon vomiting., Aspiration may cause pulmonary edema and pneumonitis.

##### *Acute inhalation toxicity*

LC50 rat: 12.9 mg/l; 1 h

OECD Test Guideline 403

Symptoms: mucosal irritations, Cough, Shortness of breath, Inhalation may lead to the formation of oedemas in the respiratory tract., Possible damages:, damage of respiratory tract

Irritating to respiratory system.

##### *Acute dermal toxicity*

LD50 rabbit: 1,530 mg/kg  
(RTECS)

---

MATERIAL SAFETY DATA SHEET  
according to the Hazard Communication Standard (29 CFR 1910.1200)

Product number 818737  
Product name Chlorotrimethylsilane for synthesis

---

Version 1.0

*Skin irritation*

rabbit

Result: Causes burns.  
(IUCLID)

Causes severe burns.

*Eye irritation*

rabbit

Result: Causes burns.  
(IUCLID)

Causes serious eye damage.

Risk of blindness!

*Genotoxicity in vivo*

Mutagenicity (mammal cell test):

Result: negative  
(IUCLID)

*Genotoxicity in vitro*

Ames test

Result: negative

Method: OECD Test Guideline 471

*Specific target organ systemic toxicity - single exposure*

May cause respiratory irritation.

*Specific target organ systemic toxicity - repeated exposure*

The substance or mixture is not classified as specific target organ toxicant, repeated exposure.

*Aspiration hazard*

Regarding the available data the classification criteria are not fulfilled.

**Carcinogenicity**

IARC	No ingredient of this product present at levels greater than or equal to 0.1% is identified as probable, possible or confirmed human carcinogen by IARC.
OSHA	No ingredient of this product present at levels greater than or equal to 0.1% is identified as a carcinogen or potential carcinogen by OSHA.
NTP	No ingredient of this product present at levels greater than or equal to 0.1% is identified as a known or anticipated carcinogen by NTP.
ACGIH	No ingredient of this product present at levels greater than or equal to 0.1% is identified as a carcinogen or potential carcinogen by ACGIH.

**Further information**

Decomposition of the substance with tissue moisture.

Other dangerous properties can not be excluded.

Handle in accordance with good industrial hygiene and safety practice.

# MATERIAL SAFETY DATA SHEET

according to the Hazard Communication Standard (29 CFR 1910.1200)

Product number 818737  
Product name Chlorotrimethylsilane for synthesis

Version 1.0

## SECTION 12. Ecological information

### Ecotoxicity

#### *Toxicity to fish*

LC0 Danio rerio (zebra fish):  $\geq 1,000$  mg/l; 96 h (External MSDS)

#### *Toxicity to bacteria*

EC50 activated sludge:  $> 10,000$  mg/l (External MSDS)

### Persistence and degradability

No information available.

### Bioaccumulative potential

No information available.

### Mobility in soil

No information available.

### Other adverse effects

#### *Additional ecological information*

Discharge into the environment must be avoided.

## SECTION 13. Disposal considerations

The information presented only applies to the material as supplied. The identification based on characteristic(s) or listing may not apply if the material has been used or otherwise contaminated. It is the responsibility of the waste generator to determine the toxicity and physical properties of the material generated to determine the proper waste identification and disposal methods in compliance with applicable regulations. Disposal should be in accordance with applicable regional, national and local laws and regulations.

## SECTION 14. Transport information

### Land transport (DOT)

UN number	UN 1298
Proper shipping name	TRIMETHYLCHLOROSILANE
Class	3 ( 8)
Packing group	II
Environmentally hazardous	--

### Air transport (IATA)

UN number	UN 1298
Proper shipping name	TRIMETHYLCHLOROSILANE
Class	3 ( 8)
Packing group	II
Environmentally hazardous	--
Special precautions for user	yes
IATA ( Passenger)	Not permitted for transport

### Sea transport (IMDG)



MATERIAL SAFETY DATA SHEET  
according to the Hazard Communication Standard (29 CFR 1910.1200)

Product number 818737  
Product name Chlorotrimethylsilane for synthesis

---

Version 1.0

UN number	UN 1298
Proper shipping name	TRIMETHYLCHLOROSILANE
Class	3 ( 8)
Packing group	II
Environmentally hazardous	--
Special precautions for user	yes
EmS	F-E S-C

---

## SECTION 15. Regulatory information

### United States of America

#### OSHA Hazards

Flammable Liquid  
Harmful by skin absorption.  
Corrosive to skin  
Corrosive to eyes  
Respiratory irritant

This information is based on 29 CFR 1910.1200 criteria prior to adoption of the GHS, and may deviate from the GHS information on the label and in section 2.

#### SARA 311/312 Hazards

Fire Hazard  
Acute Health Hazard

#### SARA 313

SARA 313: This material does not contain any chemical components with known CAS numbers that exceed the threshold (De Minimis) reporting levels established by SARA Title III, Section 313.

#### Clean Water Act

This product does not contain any Hazardous Substances listed under the U.S. CleanWater Act, Section 311, Table 116.4A.

This product does not contain any Hazardous Chemicals listed under the U.S. CleanWater Act, Section 311, Table 117.3.

### US State Regulations

#### Massachusetts Right To Know

*Ingredients*  
chlorotrimethylsilane

#### Pennsylvania Right To Know

*Ingredients*  
chlorotrimethylsilane

#### New Jersey Right To Know

*Ingredients*  
chlorotrimethylsilane

#### California Prop 65 Components

This product does not contain any chemicals known to the State of California to cause cancer, birth, or any other reproductive defects.

#### Notification status

# MATERIAL SAFETY DATA SHEET

according to the Hazard Communication Standard (29 CFR 1910.1200)

Product number	818737	Version 1.0
Product name	Chlorotrimethylsilane for synthesis	

---

TSCA: On TSCA Inventory

DSL: All components of this product are on the Canadian DSL.

---

## SECTION 16. Other information

### Training advice

Provide adequate information, instruction and training for operators.

### Full text of H-Statements referred to under sections 2 and 3.

H225	Highly flammable liquid and vapor.
H314	Causes severe skin burns and eye damage.
H335	May cause respiratory irritation.

### Key or legend to abbreviations and acronyms used in the safety data sheet

Used abbreviations and acronyms can be looked up at [www.wikipedia.org](http://www.wikipedia.org).

Date of issue: 12/20/2012

---

The information contained herein is based on the present state of our knowledge. It characterizes the product with regard to appropriate safety precautions. It does not represent a warranty of any product properties and we assume no liability for any loss or injury which may result from the use of this information. Users should conduct their own investigations to determine the suitability of the information.

*All rights reserved. Millipore and the "M" Mark are registered trademarks of Merck KGaA, Darmstadt, Germany.*