



# MATERIAL SAFETY DATA SHEET

according to the (US) Hazard Communication Standard (29 CFR 1910.1200)

Date of issue: 02/04/2013

Version 1.0

## SECTION 1. Identification

### Product identifier

Product number	808146
Product name	1-Tetradecanol for synthesis

### Relevant identified uses of the substance or mixture and uses advised against

Identified uses	Chemical for synthesis
-----------------	------------------------

### Details of the supplier of the safety data sheet

Company	EMD Millipore Corporation   290 Concord Road, Billerica, MA 01821, United States of America   SDS Phone Support: +1-978-715-1335   General Inquiries: +1-978-715-4321   Monday to Friday, 9:00 AM to 4:00 PM Eastern Time (GMT-5)  e-mail: mm_sds@merckgroup.com
---------	--

Emergency telephone	800-424-9300 CHEMTREC (USA) +1-703-527-3887 CHEMTREC (International) 24 Hours/day; 7 Days/week
---------------------	--

## SECTION 2. Hazards identification

### GHS-Labeling

Not a dangerous substance according to GHS.

### OSHA Hazards

While this material is not considered hazardous by the OSHA Hazard Communication Standard (29 CFR 1910.1200), this MSDS contains valuable information critical to the safe handling and proper use of the product. This MSDS should be retained and available for employees and other users of this product.

### Other hazards

None known.

## SECTION 3. Composition/information on ingredients

Formula	C <sub>14</sub> H <sub>30</sub> O (Hill)
CAS-No.	112-72-1
Molar mass	214.39 g/mol

# MATERIAL SAFETY DATA SHEET

according to the (US) Hazard Communication Standard (29 CFR 1910.1200)

Product number

808146

Version 1.0

Product name

1-Tetradecanol for synthesis

Remarks

No hazardous ingredients according to the OSHA Hazard Communication Standard 29 CFR 1910.1200.

---

## SECTION 4. First aid measures

### Description of first-aid measures

#### *Inhalation*

After inhalation: fresh air.

#### *Skin contact*

After skin contact: wash off with plenty of water. Remove contaminated clothing.

#### *Eye contact*

After eye contact: rinse out with plenty of water with the eyelid held wide open. Call in ophthalmologist if necessary.

#### *Ingestion*

After swallowing: make victim drink water (two glasses at most). Consult doctor if feeling unwell.

Never give anything by mouth to an unconscious person.

### Most important symptoms and effects, both acute and delayed

The following applies to aliphatic alcohols in general: effect when product is not handled and used properly: mucosal irritations; after absorption of large quantities: narcosis.

### Indication of any immediate medical attention and special treatment needed

No information available.

---

## SECTION 5. Fire-fighting measures

### Extinguishing media

#### *Suitable extinguishing media*

Water, Carbon dioxide (CO<sub>2</sub>), Foam, Dry powder

#### *Unsuitable extinguishing media*

For this substance/mixture no limitations of extinguishing agents are given.

### Special hazards arising from the substance or mixture

Combustible material, Vapors are heavier than air and may spread along floors.

Forms explosive mixtures with air on intense heating.

Development of hazardous combustion gases or vapors possible in the event of fire.

### Advice for firefighters

#### *Special protective equipment for fire-fighters*

In the event of fire, wear self-contained breathing apparatus.

---

## SECTION 6. Accidental release measures

### Personal precautions, protective equipment and emergency procedures

Advice for non-emergency personnel: Avoid generation of dusts; do not inhale dusts. Evacuate the danger area, observe emergency procedures, consult an expert.

Advice for emergency responders: Protective equipment see section 8.

# MATERIAL SAFETY DATA SHEET

according to the (US) Hazard Communication Standard (29 CFR 1910.1200)

Product number 808146  
Product name 1-Tetradecanol for synthesis

---

Version 1.0

## Environmental precautions

No special precautionary measures necessary.

## Methods and materials for containment and cleaning up

Observe possible material restrictions (see sections 7 and 10).

Take up dry. Dispose of properly. Clean up affected area. Avoid generation of dusts.

---

## SECTION 7. Handling and storage

### Precautions for safe handling

Observe label precautions.

### Conditions for safe storage, including any incompatibilities

Tightly closed. Dry.

Store at +15°C to +25°C (+59°F to +77°F).

---

## SECTION 8. Exposure controls/personal protection

### Exposure limit(s)

Contains no substances with occupational exposure limit values.

### Engineering measures

Technical measures and appropriate working operations should be given priority over the use of personal protective equipment.

### Individual protection measures

Protective clothing should be selected specifically for the workplace, depending on concentration and quantity of the hazardous substances handled. The chemical resistance of the protective equipment should be inquired at the respective supplier.

### *Hygiene measures*

Change contaminated clothing. Wash hands after working with substance.

### *Eye/face protection*

Safety glasses

### *Hand protection*

Chemical-resistant, impervious gloves complying with an approved standard should be worn at all times when handling chemical products if a risk assessment indicates this is necessary.

### *Respiratory protection*

required when dusts are generated.

Use a properly fitted, air-purifying or air-fed respirator complying with an approved standard if a risk assessment indicates this is necessary. Respirator selection must be based on known or anticipated exposure levels, the hazards of the product and the safe working limits of the selected respirator.

---

## SECTION 9. Physical and chemical properties

Physical state	scales
Color	white
Odor	fatty odor

# MATERIAL SAFETY DATA SHEET

according to the (US) Hazard Communication Standard (29 CFR 1910.1200)

Version 1.0

Product number 808146  
Product name 1-Tetradecanol for synthesis

Odor Threshold	No information available.
pH	No information available.
Melting point	36 - 39 °C
Boiling point/boiling range	531 - 550 °F ( 277 - 288 °C) at 1,013 hPa
Flash point	293 °F ( 145 °C) Method: DIN 51758
Evaporation rate	No information available.
Flammability (solid, gas)	No information available.
Lower explosion limit	No information available.
Upper explosion limit	No information available.
Vapor pressure	< 1 hPa at 68 °F ( 20 °C)
Relative vapor density	No information available.
Relative density	0.8 g/cm <sup>3</sup> at 140 °F ( 60 °C)
Water solubility	at 68 °F ( 20 °C) practically insoluble
Partition coefficient: n-octanol/water	log Pow: 6.03 (experimental) (Lit.) A remarkable bioaccumulation potential is expected (log Po/w >3).
Autoignition temperature	No information available.
Decomposition temperature	No information available.
Viscosity, dynamic	6.4 mPa.s at 140 °F ( 60 °C)
Explosive properties	No information available.
Ignition temperature	464 °F ( 240 °C) Method: DIN 51794

## SECTION 10. Stability and reactivity

### Reactivity

## MATERIAL SAFETY DATA SHEET

according to the (US) Hazard Communication Standard (29 CFR 1910.1200)

Product number 808146  
Product name 1-Tetradecanol for synthesis

---

Version 1.0

Forms explosive mixtures with air on intense heating.

The following applies in general to flammable organic substances and mixtures: in correspondingly fine distribution, when whirled up a dust explosion potential may generally be assumed.

### Chemical stability

The product is chemically stable under standard ambient conditions (room temperature) .

### Possibility of hazardous reactions

Violent reactions possible with:

strong alkalis, Strong oxidizing agents, Acid chlorides, Acid anhydrides

### Conditions to avoid

Strong heating.

### Incompatible materials

no information available

### Hazardous decomposition products

no information available

---

## SECTION 11. Toxicological information

### Information on toxicological effects

#### *Likely route of exposure*

Eye contact, Skin contact, Ingestion

#### *Acute oral toxicity*

LD50 rat: > 2,000 mg/kg

OECD Test Guideline 401

#### *Acute dermal toxicity*

LD50 rabbit: > 2,000 mg/kg

OECD Test Guideline 402

#### *Skin irritation*

rabbit

Result: No irritation

OECD Test Guideline 404

#### *Eye irritation*

rabbit

Result: slight irritation

OECD Test Guideline 405

#### *Sensitization*

Patch test: human

Result: negative

(External MSDS)

#### *Genotoxicity in vitro*

Ames test

Result: negative

Method: OECD Test Guideline 471

## MATERIAL SAFETY DATA SHEET

according to the (US) Hazard Communication Standard (29 CFR 1910.1200)

Product number 808146  
Product name 1-Tetradecanol for synthesis

---

Version 1.0

### *Specific target organ systemic toxicity - single exposure*

The substance or mixture is not classified as specific target organ toxicant, single exposure.

### *Specific target organ systemic toxicity - repeated exposure*

The substance or mixture is not classified as specific target organ toxicant, repeated exposure.

### *Aspiration hazard*

Regarding the available data the classification criteria are not fulfilled.

## **Carcinogenicity**

IARC	No ingredient of this product present at levels greater than or equal to 0.1% is identified as probable, possible or confirmed human carcinogen by IARC.
OSHA	No ingredient of this product present at levels greater than or equal to 0.1% is identified as a carcinogen or potential carcinogen by OSHA.
NTP	No ingredient of this product present at levels greater than or equal to 0.1% is identified as a known or anticipated carcinogen by NTP.
ACGIH	No ingredient of this product present at levels greater than or equal to 0.1% is identified as a carcinogen or potential carcinogen by ACGIH.

## **Further information**

Hazardous properties cannot be excluded but are unlikely when the product is handled appropriately.

### Other information

The following applies to aliphatic alcohols in general: effect when product is not handled and used properly: mucosal irritations; after absorption of large quantities: narcosis.

### Further data:

Handle in accordance with good industrial hygiene and safety practice.

---

## **SECTION 12. Ecological information**

### **Ecotoxicity**

#### *Toxicity to fish*

LC50 Danio rerio (zebra fish): > 10,000 mg/l; 96 h (above the solubility limit in the test medium) (IUCLID)

#### *Toxicity to bacteria*

EC50 Pseudomonas putida: > 10,000 mg/l; 30 min

OECD Test Guideline 209 (above the solubility limit in the test medium)

### **Persistence and degradability**

#### *Biodegradability*

> 60 %; 28 d

OECD Test Guideline 301D

Readily biodegradable.

### **Bioaccumulative potential**

## MATERIAL SAFETY DATA SHEET

according to the (US) Hazard Communication Standard (29 CFR 1910.1200)

Product number 808146  
Product name 1-Tetradecanol for synthesis

---

Version 1.0

*Partition coefficient: n-octanol/water*

log Pow: 6.03

(experimental)

(Lit.) A remarkable bioaccumulation potential is expected (log Po/w >3).

### **Mobility in soil**

No information available.

### **Other adverse effects**

*Additional ecological information*

No ecological problems are to be expected when the product is handled and used with due care and attention.

---

## **SECTION 13. Disposal considerations**

The information presented only applies to the material as supplied. The identification based on characteristic(s) or listing may not apply if the material has been used or otherwise contaminated. It is the responsibility of the waste generator to determine the toxicity and physical properties of the material generated to determine the proper waste identification and disposal methods in compliance with applicable regulations. Disposal should be in accordance with applicable regional, national and local laws and regulations.

---

## **SECTION 14. Transport information**

### **Land transport (DOT)**

Not classified as dangerous in the meaning of transport regulations.

### **Air transport (IATA)**

Not classified as dangerous in the meaning of transport regulations.

### **Sea transport (IMDG)**

Not classified as dangerous in the meaning of transport regulations.

---

## **SECTION 15. Regulatory information**

### **United States of America**

#### **OSHA Hazards**

No OSHA Hazards

This information is based on 29 CFR 1910.1200 criteria prior to adoption of the GHS, and may deviate from the GHS information on the label and in section 2.

#### **SARA 311/312 Hazards**

No SARA Hazards

#### **SARA 313**

SARA 313: This material does not contain any chemical components with known CAS numbers that exceed the threshold (De Minimis) reporting levels established by SARA Title III, Section 313.

## MATERIAL SAFETY DATA SHEET

according to the (US) Hazard Communication Standard (29 CFR 1910.1200)

Product number 808146  
Product name 1-Tetradecanol for synthesis

---

Version 1.0

### SARA 302

SARA 302: No chemicals in this material are subject to the reporting requirements of SARA Title III, Section 302.

### Clean Water Act

This product does not contain any Hazardous Substances listed under the U.S. CleanWater Act, Section 311, Table 116.4A.

This product does not contain any Hazardous Chemicals listed under the U.S. CleanWater Act, Section 311, Table 117.3.

### Massachusetts Right To Know

Remarks

No components are subject to the Massachusetts Right to Know Act.

### Pennsylvania Right To Know

*Ingredients*

1-tetradecanol, n-

### New Jersey Right To Know

*Ingredients*

1-tetradecanol, n-

1-Dodecanol

### Notification status

TSCA: On TSCA Inventory

DSL: All components of this product are on the Canadian DSL.

---

## SECTION 16. Other information

### Training advice

Provide adequate information, instruction and training for operators.

### Key or legend to abbreviations and acronyms used in the safety data sheet

Used abbreviations and acronyms can be looked up at [www.wikipedia.org](http://www.wikipedia.org).

Date of issue:02/04/2013

---

The information contained herein is based on the present state of our knowledge. It characterizes the product with regard to appropriate safety precautions. It does not represent a warranty of any product properties and we assume no liability for any loss or injury which may result from the use of this information. Users should conduct their own investigations to determine the suitability of the information.

*All rights reserved. Millipore and the "M" Mark are registered trademarks of Merck KGaA, Darmstadt, Germany.*