

## Datasheet

### GM polyclonal antibody

**Catalog Number:** PAB8280

**Regulation Status:** For research use only (RUO)

**Product Description:** Sheep polyclonal antibody raised against recombinant GM-CSF.

**Immunogen:** Recombinant protein corresponding to human GM-CSF.

**Host:** Sheep

**Reactivity:** Human

**Applications:** EIA, WB

(See our web site product page for detailed applications information)

**Protocols:** See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Form:** Liquid

**Recommend Usage:** The optimal working dilution should be determined by the end user.

**Storage Buffer:** In 20 mM phosphate buffer, 150 mM NaCl, pH 7.2 (0.09% sodium azide)

**Storage Instruction:** Store at -20°C.

Aliquot to avoid repeated freezing and thawing.