

Datasheet

TMOD4 MaxPab mouse polyclonal antibody (B01)

Catalog Number: H00029765-B01

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length human TMOD4 protein.

Immunogen: TMOD4 (AAH17810.1, 1 a.a. ~ 345 a.a) full-length human protein.

Sequence:

MSSYQKELEKYRDIDEDEILRTLSPEELEQLDCEIQEM
DPENMLLPAGLRQRDQTKKSPTGPLDREALLQYLEQQ
ALEVKERDDLVPFTGEKKGKPYIQPKREIPAAEQITLPE
ELEEALAHATDAEMCDIAAILDMYTLMSNKQYYDALCS
GEICNTEGISSVVQPDKYKVPDEPPNPTNIEEILKRVR
SNDKELEEVLNNIQDIPMLSELCEAMKANTYVRSF
SLVATRSGDPIANAVADMLRENRLQSLNIESNFISST
GLMAVLKAVRENATLTCLRVDNQRQWPGDAVEMEMA
TVLEQRPSIVRFGYHFTQQGPRARAAQAMTRNNELRR
QQKKR

Host: Mouse

Reactivity: Human

Applications: WB-Ce, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: No additive

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 29765

Gene Symbol: TMOD4

Gene Alias: DKFZp779I0852, SK-TMOD