

Datasheet

SNX8 monoclonal antibody (M08), clone 3B6

Entrez GeneID: 29886**Gene Symbol:** SNX8**Gene Alias:** Mvp1**Catalog Number:** H00029886-M08**Regulation Status:** For research use only (RUO)**Product Description:** Mouse monoclonal antibody
raised against a full-length recombinant SNX8.**Clone Name:** 3B6**Immunogen:** SNX8 (AAH21565, 1 a.a. ~ 465 a.a)
full-length recombinant protein with GST tag. MW of the
GST tag alone is 26 KDa.**Sequence:**

MTGRAMDPLPAAAVGAAAEAEDEEADPPASDLPTP
QAIEPQAI VQQVPAPSRMQMPQGNPLLSHTLQELLA
RDTVQVELIPEKKGLFLKHVEYEVSSQRFKSSVYRRYN
DFVVFQEMLLHKFPYRMVPALPPKRMLGADREFIEAR
RRALKRFVNLVARHPLFSEDVVLKFLSFSGSDVQNKL
KESAQCVGDEFLNCKLATRAKDFLPADIQAQFAISRELI
RNIYNSFHKLRDRAERIASRAIDNAADLLIFGKELSAIGS
DTTPLPSWAALNSSTWGS LKQALKGLSVEFALLADKA
AQQGKQEENDVVEKLNFLDLLQSYKDL CERHEKGVL
HKHQALHKYSLMKRQMMSATAQNREPESVEQLESR
IVEQENAIQTMELRNYFSLYCLHQETQLIHVYLP L TSHIL
RAFVNSQIQGHKEMSKVWNDLRPKLSCLFAGPHSTLT
PPCSPPEDGLCPH

Host: Mouse**Reactivity:** Human,Rat**Applications:** ELISA, S-ELISA, WB-Ce(See our web site product page for detailed applications
information)**Protocols:** See our web site at<http://www.abnova.com/support/protocols.asp> or product
page for detailed protocols**Isotype:** IgG2b Kappa**Storage Buffer:** In 1x PBS, pH 7.4**Storage Instruction:** Store at -20°C or lower. Aliquot to
avoid repeated freezing and thawing.