

Datasheet

SYCP3 purified MaxPab rabbit polyclonal antibody (D01P)

Catalog Number: H00050511-D01P

Regulation Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against a full-length human SYCP3 protein.

Immunogen: SYCP3 (NP_710161.1, 1 a.a. ~ 236 a.a) full-length human protein.

Sequence:

MVSSGKKYSRKSGKPSVEDQFTRAYDFETEDKKDLS
GSEEDVIEGKTAVIEKRRKKRSSAGVVEDMGGEVQNM
LEGVGVDINKALLAKRKRLEMYTKASLKTSNQKIEHVV
KTQQDQRQKLNQEYSQQFLTLFQQWDLDMQKAEQ
EEKILNMFRQQKILQQSRIVQSQRKTIKQLYEQFIKS
MEELEKNHDNLLTGAQNEFKKEMAMLQKKIMMETQQ
QEIASVRKSLQSMLF

Host: Rabbit

Reactivity: Human

Applications: WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 50511

Gene Symbol: SYCP3

Gene Alias: COR1, MGC71888, SCP3