

Datasheet

WDR42A monoclonal antibody (M01A), clone 1A12

Catalog Number: H00050717-M01A

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against a full length recombinant WDR42A.

Clone Name: 1A12

Immunogen: WDR42A (AAH13107, 1 a.a. ~ 273 a.a)
full-length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.

Sequence:

MSSKGSSTDGRTDLANGSLSSSPEEMSGAEEGRETS
SGIEVEASDLSLSLTGDDGGPNRTSTESRGTESSG
EDKSDSMEDTGHYSINDENRVHDERSEEEEEEEEEEE
EEQPRRRRVQRKRANRDQDSSDDERALEDWVSSETSA
LPRPRWQALPALRERELGSSARFVYEACGARVVFQR
FRLQHGLEGHTGCVNTLHFNQRGTWLASGSDDLKVV
VWDWVRRQPVLDFESGHKSNVFQVRQGSIIATERIRH
ELEKSELQHGLGADSELL

Host: Mouse

Reactivity: Human

Applications: ELISA, WB-Ce, WB-Re

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Isotype: IgM kappa

Storage Buffer: In ascites fluid

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 50717

Gene Symbol: WDR42A

Gene Alias: DKFZp781G1096, FLJ35857, H326,

MGC117276, MGC118891, MGC99640

Gene Summary: This gene encodes a WD repeat-containing protein that interacts with the Cul4-Ddb1 E3 ligase macromolecular complex. Multiple alternatively spliced transcript variants have been found for this gene. [provided by RefSeq]