

Datasheet

SPOCK3 polyclonal antibody (A01)

Catalog Number: H00050859-A01

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a partial recombinant SPOCK3.

Immunogen: SPOCK3 (NP_058646, 132 a.a. ~ 230 a.a) partial recombinant protein with GST tag.

Sequence:

RGPISTCKQCPVVYPSPVCGSDGHTYSFQCKLEYQA
CVLGKQISVKCEGHCPSPDKPTSTSRNVKRACSDLE
FREVANRLRDWFKALHESGSQNKKT

Host: Mouse

Reactivity: Human

Applications: ELISA, WB-Ce, WB-Re

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: 50 % glycerol

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 50859

Gene Symbol: SPOCK3

Gene Alias: HSAJ1454, TES-3

Gene Summary: Proteoglycans, which consist of a core protein and covalently linked glycosaminoglycans, are components of the extracellular matrix. SPOCK3 encodes a member of a novel Ca(2+)-binding proteoglycan family.[supplied by OMIM]