

## Datasheet

### GLRX2 polyclonal antibody (A01)

**Catalog Number:** H00051022-A01

**Regulation Status:** For research use only (RUO)

**Product Description:** Mouse polyclonal antibody raised against a full-length recombinant GLRX2.

**Immunogen:** GLRX2 (AAH28113, 1 a.a. ~ 124 a.a)  
full-length recombinant protein with GST tag.

**Sequence:**

MESNTSSSLENLATAPVNQIQETISDNCVVIFSKTSCSY  
CTMAKKLFHDMNVNYKVVELDLLEYGNQFQDALYKMT  
GERTVPRIFVNGTFIGGATDTHRLHKEGKLLPLVHQCY  
LKSKSRKEFQ

**Host:** Mouse

**Reactivity:** Human

**Applications:** ELISA, WB-Re

(See our web site product page for detailed applications information)

**Protocols:** See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Storage Buffer:** 50 % glycerol

**Storage Instruction:** Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

**Entrez GeneID:** 51022

**Gene Symbol:** GLRX2

**Gene Alias:** GRX2, bA101E13.1

**Gene Summary:** Glutaredoxins (e.g., GLRX; MIM 600443) are a family of glutathione-dependent hydrogen donors that participate in a variety of cellular redox reactions.[supplied by OMIM]