

Datasheet

BOLA1 purified MaxPab mouse polyclonal antibody (B01P)

Catalog Number: H00051027-B01P

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length human BOLA1 protein.

Immunogen: BOLA1 (NP_057158.1, 1 a.a. ~ 137 a.a) full-length human protein.

Sequence:

MLSGRLVLGLVSMAGRVCLCQGSAGSGAIGPVEAAIR
TKLEEALSPVLELRNESGGHAVPPGSETHFRVAVVS
SRFEGLSPLQRHRLVHAALAEELGGPVHALAIQARTPA
QWRENSQLDTSPCLGGNKKTLGTP

Host: Mouse

Reactivity: Human

Applications: IF, WB-Ce, WB-Tr
(See our web site product page for detailed applications information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 51027

Gene Symbol: BOLA1

Gene Alias: CGI-143, MGC75015, RP11-196G18.18