

Datasheet

ZNF593 purified MaxPab mouse polyclonal antibody (B01P)

Catalog Number: H00051042-B01P

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length human ZNF593 protein.

Immunogen: ZNF593 (AAH02580, 1 a.a. ~ 116 a.a)
full-length human protein.

Sequence:

MKAKRRRRLDEIHRELRPQGSARPQDPNAEFDPDL
PGGGLHRCLACARYFIDSANLKTFRSKDHKKRLKQL
SVEPYSQEEAERAAGMGSYVPPRR LAVPTEVSTEVPE
MDTST

Host: Mouse

Reactivity: Human

Applications: IF, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 51042

Gene Symbol: ZNF593

Gene Alias: ZT86