

Datasheet

PI15 purified MaxPab mouse polyclonal antibody (B01P)

thus may play a role in the central nervous system.
[provided by RefSeq]

Catalog Number: H00051050-B01P

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length human PI15 protein.

Immunogen: PI15 (NP_056970.1, 1 a.a. ~ 258 a.a)
full-length human protein.

Sequence:

MIAISAVSSALLFSLCEASTVLLNSTDSSPPTNNFTDI
EAALKAQLDSADIPKARRKRYISQNDMIAILDYHNQVR
GKVFPAAANMEYMWVDENLAKSAEAWAATCIWDHGP
SYLLRFLGQNL SVRTGRYRSILQLVKPWYDEVKDYAF
PYPQDCNPRCPMRCFGPMCTHYTQM VWATSNRIGC
AIHTCQNMNVWGSVWRRAYLVLCNYAPKGNWIGEAP
YKVGVPCCSSCPPSYGGSCDNLCPFVTSNLYWFK

Host: Mouse

Reactivity: Human

Applications: WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 51050

Gene Symbol: PI15

Gene Alias: CRISP8, DKFZp686F0366, P24TI, P25TI

Gene Summary: This gene encodes a trypsin inhibitor. The protein shares similarity to insect venom allergens, mammalian testis-specific proteins and plant pathogenesis-related proteins. It is frequently expressed in human neuroblastoma and glioblastoma cell lines, and