

Datasheet

SH3GLB1 purified MaxPab rabbit polyclonal antibody (D01P)

Catalog Number: H00051100-D01P

Regulation Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against a full-length human SH3GLB1 protein.

Immunogen: SH3GLB1 (NP_057093.1, 1 a.a. ~ 365 a.a) full-length human protein.

Sequence:

MNIMDFNVKKLAADAGTFLSRAVQFTEEKLGQAEKTE
LDAHLENLLSKAECTKIWTEKIMKQTEVLLQPNPNARIE
EFVYEKLDKAPSRINNPELLGQYMIDAGTEFGPGTAY
GNALIKCGETQKRIGTADRELIQTSALNFLTPLRNFIEG
DYKTIAKERKLLQNKRLDLDAKTRLKKAKAAETRNS
EQELRITQSEFDRQAEITRLLLEGISSTHAHHLRCLNDF
VEAQMTYYAQCYYMLDLQKQLGSFSPNYLSNNNQ
SVTPVPSVLPNAIGSSAMASTSGLVITSPSNLSDLKEC
SGSRKARVLYDYDAANSTELSLLADEVITVFSVVGMD
DWLMGERGNQKGKVPITYLELLN

Host: Rabbit

Reactivity: Human

Applications: WB-Ce, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 51100

Gene Symbol: SH3GLB1

Gene Alias: Bif-1, CGI-61, KIAA0491, dJ612B15.2