

Datasheet

TFB1M MaxPab mouse polyclonal antibody (B01)

Catalog Number: H00051106-B01

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length human TFB1M protein.

Immunogen: TFB1M (AAH17788.1, 1 a.a. ~ 346 a.a)
full-length human protein.

Sequence:

MAASGKLSTCRLPPLPTIREIILRLQAQKLSQNFLL
DLRLTDKIVRKAGNLNAYVYEVGPGPGGITSILNAD
VAELLVVEKDTRFIPGLQMLSDAAPGKLRIHGDVLT
KVEKAFSESLKRPWEDDPPNVHIIGNLPFSVSTPLIKW
LENISCRDGPVYGRQTMTLTFQKEVAERLAANTGSK
QRSRLSVMAQYLCNVRHIFTIPGQAFVPKPEVDVGVV
HFTPLIQPKIEQPFKLVEKVQNVFQFRRKYCHRGLRM
LFPEAQRLESTGRLLLEADIDPTLRPRQLSISFKSLCD
VYRKMCEDEPQLFAYNFREELKRRKSKNEEKEEDDA
ENYRL

Host: Mouse

Reactivity: Human

Applications: Det Ab, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: No additive

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 51106

Gene Symbol: TFB1M

Gene Alias: CGI-75, CGI75, mtTFB, mtTFB1

Gene Summary: The transcription of genes from mitochondrial DNA requires a mitochondrial RNA

polymerase (see POLRMT, MIM 601778) and a DNA-binding transcription factor (see TFAM, MIM 600438). Transcription factor B1 (TFB1M) is a part of this transcription complex.[supplied by OMIM]