

Datasheet

DBR1 purified MaxPab mouse polyclonal antibody (B01P)

Catalog Number: H00051163-B01P

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length human DBR1 protein.

Immunogen: DBR1 (NP_057300.2, 1 a.a. ~ 544 a.a) full-length human protein.

Sequence:

MRVAVAGCCHGELDKIYETLALAERRGPGPVDLLCC
GDFQAVRNEADLRCAVPPKYRHMQTFRYYSGEKK
APVLTFIGGNHEASNHLQELPYGGWVAPNIYYLGLAG
VVKYRGVRIGGISGIFKSHDYRKGFECPPYNSSTIRSI
YHVRNIEVYKLKQLKQPIDIFLSHDWPRSIYHYGNKKQL
LTKTSFFRQEVENNTLGSPAASELLEHLKPTYWFS AHL
HVKFAALMQHQAKDKGQTARATKFLALDKCLPHRDFL
QILEIEHDPSAPDYLEYDIEWLTILRATDDLINVTGRLW
NMPENNGLHARWDYSATEEGMKEVLEKLNHDLKVPC
NFSVTAACYDPSKPQTQMQLIHRINPQTTEFCAQLGIID
INVRLQKSKEEHHVCGEYEEQDDVESNDSGEDQSEY
NTDTSALSSINPDEIMLDEEEDEDSIVSAHSGMNTPSV
EPSDQASEFSASFSDVRILPGSMIVSSDDTV DSTIDRE
GKPGGTVESGNGEDLTKVPLKRLSDEHEPEQRKKIKR
RNQAIYAAVDDDDDDAA

Host: Mouse

Reactivity: Human

Applications: Det Ab, IF, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 51163

Gene Symbol: DBR1

Gene Alias: -

Gene Summary: The RNA lariat debranching enzyme, or DBR1, specifically hydrolyzes 2-prime-to-5-prime branched phosphodiester bonds at the branch point of excised lariat intron RNA and converts them into linear molecules.[supplied by OMIM]