

Datasheet

TUBD1 monoclonal antibody (M05), clone 3C7

Catalog Number: H00051174-M05

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody
raised against a partial recombinant TUBD1.

Clone Name: 3C7

Immunogen: TUBD1 (NP_057345.2, 354 a.a. ~ 453
a.a) partial recombinant protein with GST tag. MW of the
GST tag alone is 26 KDa.

Sequence:

QSADVEGFKDPALYTSWLKPVNAFNVWKTQRAFSKY
EKSAVLVSNSQFLVKPLDMIVGKAWNMFASKAYIHQY
TKFGIEEEDFLDSFTSLEQVVASYCNL

Host: Mouse

Reactivity: Human

Applications: ELISA, IP, S-ELISA

(See our web site product page for detailed applications
information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product
page for detailed protocols

Isotype: IgG2a Kappa

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to
avoid repeated freezing and thawing.

Entrez GeneID: 51174

Gene Symbol: TUBD1

Gene Alias: FLJ12709, TUBD