

Datasheet

RNF6 MaxPab mouse polyclonal antibody (B01)

Catalog Number: H00006049-B01

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length human RNF6 protein.

Immunogen: RNF6 (NP_005968.1, 1 a.a. ~ 685 a.a) full-length human protein.

Sequence:

MNQSRSRSDGGSEETLPQDHNHHENERRWQQLRH
REEAYYQFINELNDEYRLMRDHNLLGTPGEITSEELQ
QRLDGVKEQLASQPDRLDGTNYRDSEVPRESSHEDS
LLEWLNTFRRTGNATRSQGNGNQTWRAVSRTNPN
GEFRFSLEIHVNHNENRGFEIHGEDYTDIPLSDSNRDHT
ANRQQRSTSPVARRTSQTSVNFNGSSSNIPRTLAS
RGQNPAEGSFSTLGRLRNGIGGAAGIPRANASRTNFS
SHTNQSGGSELRQREGQRFGAHVWENGARSNTV
RNTNQRLPIRLRSTNSRSRSPQRQSGTVYHNSQR
ESRPVQQTTRRSVRRRGRTRVFLEQDRERERRGTAY
TPFSNSRLVSRITVEEGEESSRSSTAVRRHPTITLDLQ
VRRIRPGENRDRDSIANRTRSRVGLAENTVTIESNSGG
FRRTISRLERSGIRTYVSTITVPLRRISENELVEPSSVAL
RSILRQIMTGFGELSSLMEADSESELQRNGQHLPDMH
SELSNLGTDNNRSQHREGSSQDRQAQGDSTEMHGE
NETTQPHTRNSDSRGGRLRNPNNLVETGTLPIRLA
HFFLLNESDDDDRIRGLTKEQIDNLSTRHYEHNSIDSEL
GKICSVICISDYVTGNKLRQLPCMHEFHIHCIDRWLSEN
CTCPICRQPVLGSNIANNG

Host: Mouse

Reactivity: Human

Applications: WB-Ce, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: No additive

Storage Instruction: Store at -20°C or lower. Aliquot to

avoid repeated freezing and thawing.

Entrez GeneID: 6049

Gene Symbol: RNF6

Gene Alias: DKFZp686P0776

Gene Summary: The protein encoded by this gene contains a RING-H2 finger motif. Deletions and mutations in this gene were detected in esophageal squamous cell carcinoma (ESCC), suggesting that this protein may be a potential tumor suppressor. Studies of the mouse counterpart suggested a role of this protein in the transcription regulation that controls germinal differentiation. Multiple alternatively spliced transcript variants encoding the same protein are observed. [provided by RefSeq]