

Datasheet

HDHD1A monoclonal antibody (M05), clone 4F12

Catalog Number: H00008226-M05

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody
raised against a partial recombinant HDHD1A.

Clone Name: 4F12

Immunogen: HDHD1A (NP_036212.2, 107 a.a. ~ 212
a.a) partial recombinant protein with GST tag. MW of the
GST tag alone is 26 KDa.

Sequence:

DMKTSRHKEFFSLFSHIVLGDDPEVQHKGKPDPIFLAC
AKRFSPPPAMEKCLVFEDAPNGVEAALAAGMQVVMV
PDGNLSRDLTTKATLVLSLQDFQPELFGPLPS

Host: Mouse

Reactivity: Human

Applications: ELISA, S-ELISA, WB-Re, WB-Tr
(See our web site product page for detailed applications
information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product
page for detailed protocols

Isotype: IgG2a Kappa

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to
avoid repeated freezing and thawing.

Entrez GeneID: 8226

Gene Symbol: HDHD1A

Gene Alias: DXF68S1E, FAM16AX, GS1