

Datasheet

CSDA purified MaxPab mouse polyclonal antibody (B01P)

Catalog Number: H00008531-B01P

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length human CSDA protein.

Immunogen: CSDA (AAH15564.1, 1 a.a. ~ 372 a.a) full-length human protein.

Sequence:

MSEAGEATTTTTTLPQAPTEAAAAAPQDPAPKSPVG
SGAPQAAAPAPAAHVAGNPGGDAAPAATGTAAAASL
ATAAGSEDAEKKVLATKVLGTVKWFNVRNGYGFNIRN
DTKEDVFVHQTAIKKNNPRKYLRVGDGETVEFDVVE
GEKGAEAAANTGPDGVPVEGSRYAADRRRRYRRGY
GRRRGPPRNYAGEEEEEEGSGSSEGFDPATDRQFSG
ARNQLRRPQYRPQYRQRRFPYHVGQTFDRRSRVLP
HPNRIQAGEIGEMKDGVPGEAQLQGPVHRNPTYRPR
YRSRGPPRPRPAPAVGEAEDKENQQATSGPNQPSVR
RGYRRPYNYRRRPRPPNAPSQDGKEAKAGEAPTENP
APPTQQSSAE

Host: Mouse

Reactivity: Human

Applications: Det Ab, IF, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 8531

Gene Symbol: CSDA

Gene Alias: CSDA1, DBPA, ZONAB