

Datasheet

BLZF1 purified MaxPab mouse polyclonal antibody (B01P)

Catalog Number: H00008548-B01P

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length human BLZF1 protein.

Immunogen: BLZF1 (NP_003657.1, 1 a.a. ~ 400 a.a) full-length human protein.

Sequence:

MTTKNLETKVTVTSSPIRGAGDGMETEEPPKSVEVTS
GVQSRKHHSLLQSPWKKAVPSESPGVLQLGKMLTEKA
MEVKAVRILVPKAAITHDIPNKNKVKSLGHHKGEFLG
QSEGVIEPNKELSEVKNVLEKLKNSERRLLQDKEGLSN
QLRVQTEVNRELKLLVASVGDDLQYHFERLAREKNQ
LILENEALGRNTAQLSEQLERMSIQCDVWRSKFLASRV
MADELTNSRAALQRQNRDAHGAIQDLLSEREQFRQE
MIATQKLLLELLVSLQWGREQTYSPSVQPHSTAELALT
NHKLAKAVNSHLLGNVGINNQKKIPSTVEFCSTPAEKM
AETVLRILDPVTCKESSPDNPFESSPTLLATKKNIGR
FHPYTRYENITFNCCNHCRGELIAL

Host: Mouse

Reactivity: Human

Applications: IF, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 8548

Gene Symbol: BLZF1

Gene Alias: GOLGIN-45, JEM-1, JEM-1s, JEM1, MGC22497