

Datasheet

SCRT1 polyclonal antibody

Catalog Number: PAB24402

Regulation Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against recombinant SCRT1.

Immunogen: Recombinant protein corresponding to amino acids of human SCRT1.

Sequence:

LHDKGYLSDYVGPSSVYDGDAAALLKGPSPEPMYAA
AVRGELGP

Host: Rabbit

Reactivity: Human

Applications: IHC-P, WB

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Purification: Antigen affinity purification

Isotype: IgG

Recommend Usage: Immunohistochemistry (1:20-1:50)

Western Blot (1:250-1:500)

The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS, pH 7.5 (40% glycerol, 0.02% sodium azide)

Storage Instruction: Store at 4°C. For long term storage store at -20°C.

Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 83482

Gene Symbol: SCRT1

Gene Alias: DKFZp547F072, SCRT

Gene Summary: This gene is a member of the Snail family of C2H2-type zinc finger transcription factors. It codes for a neural-specific transcriptional repressor that binds to E-box motifs. The protein may promote neural differentiation and may be involved in cancers with neuroendocrine features. [provided by RefSeq]