



SCI-FAB

Quality stainless steel furnishings for research environments

Acid Cabinet Base Unit



Overall Dimensions			
Product Number	W	D	H
BAC242235	24"	22"	34 3/4"
• Flush front panel above			

Overall Dimensions			
Product Number	W	D	H
BAC302235	30"	22"	34 3/4"
BAC362235	36"	22"	34 3/4"
BAC422235	42"	22"	34 3/4"
BAC482235	48"	22"	34 3/4"
• Flush front panel above			

Common Features

- SEFA section 8 tested and approved
- Spill containment trays
- "ACID" cabinet labeling
- ADA height available
- Adjustable shelf
- Louvered doors
- Lock options available



Phone: 1-800-346-3246

Fax: (303) 789-1741

E-mail: dhall@kloppenberg.com

www.kloppenberg.com