
Product Specifications

VWR Part Number: IR1358
ImmunoReagents, Inc. Part Number: GtxOt-070-DALP

Product: GOAT ANTI-BIO ALKPPOS 1MG

Concentration: 1.50 mg/ml (E 1% at 280 nm = 13.0)

Amount: 1.0 mg

Conjugate: Alkaline Phosphatase

Form: Clear, colorless liquid

Purification: Affinity purified using solid phase biotin

Purity: > 95% based on SDS-PAGE

Host: Goat

Immunogen: BSA-Biotin, KLH-Biotin

Buffer: 30 mM Triethanolamine, pH 7.2, 5 mM Magnesium Chloride, 0.1 mM Zinc Chloride, 1 % (w/v) BSA, Protease/IgG free

Storage: Store undiluted liquid at 2-8 °C. For storage at -20 °C, dilute with an equal volume of glycerol to prevent loss of enzymatic activity. Prepare working dilutions prior to use and then discard.

Shelf Life: 1 year from date of receipt. Prepare working dilution prior to use and then discard.

Specificity: This antibody reacts with Biotin

Cross Reactivity: This antibody does not react with free BSA or free KLH

Country of Origin: Goat serum was obtained from healthy animals of US origin and under the care of a registered veterinarian.

Applications: This conjugate is suitable for all immunoassay applications. The optimal working dilution should be determined by the investigator. Suggested starting dilution(?): · 1:500-1:2,000 for Immunohistochemistry · 1:50,000-1:5,000 for ELISA/Western blot

Disclaimer: For in vitro Laboratory Use Only. Not for diagnostic or therapeutic use. Not for human or animal consumption. Suggested applications of our products are not recommendations to use our products in violation of any patent or as a license under any patent of ImmunoReagents, Inc.

Product may not be resold or modified for resale without prior written approval of ImmunoReagents, Inc.