
Product Specifications

VWR Part Number: IR2733
ImmunoReagents, Inc. Part Number: MAB-100-IL17ABIO

Product: MONOCLONAL ANTI-IL-17A

Concentration: 0.5 mg/mL solution

Amount: 100 µg

Conjugate: Biotin (NHS-LC)

Form: 0.5 mg/mL solution

Purification: Protein A

Purity: > 95% based on SDS-PAGE

Host: Mouse

Immunogen: Interleukin-17A

Buffer: 0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2 with no stabilizers or preservatives

Storage: Store at 2-8°C for use within 1-3 weeks or -20°C for long term storage, avoid repeat freeze thaws.

Shelf Life: 1 year from date of receipt. Prepare working dilution prior to use and then discard.

Specificity: Interleukin-17A

Cross Reactivity: Interleukin-17A (IL-17A), also known as CTLA-8, is a pro-inflammatory cytokine member of a six-species family of proteins (IL-17A-17F). IL-17A is secreted mainly by activated CD4+ and CD8+ T lymphocytes and acts through its receptor, IL-17R, to induce

Country of Origin: US Origin

Applications: WB

Disclaimer: For in vitro Laboratory Use Only. Not for diagnostic or therapeutic use. Not for human or animal consumption. Suggested applications of our products are not recommendations to use our products in violation of any patent or as a license under any patent of ImmunoReagents, Inc.

Product may not be resold or modified for resale without prior written approval of ImmunoReagents, Inc.