

Product Specifications

VWR Part Number: IR2743
ImmunoReagents, Inc. Part Number: PAR-100-IL4

Product: RABBIT ANTI-IL-4

Concentration: N/A

Amount: 100 µg

Conjugate: Unconjugated

Form: Lyophilized

Purification: Protein A

Purity: > 90% based on SDS-PAGE

Host: Rabbit

Immunogen: Human IL-4

Buffer: 0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2 with no stabilizers or preservatives

Storage: Store at 2-8°C after restoration with 0.1mL deionized water or -20°C for long term storage, avoid repeat freeze thaws.

Shelf Life: 1 year from date of receipt. Prepare working dilution prior to use and then discard.

Specificity: Human IL-4

Cross Reactivity: Rabbit anti-human IL-4 recognizes native and recombinant forms of human IL-4 in tissue culture supernatant and other body fluids. A control of similar fluids known to contain no human IL-4 should be used as a control.

Country of Origin: US Origin

Applications: For use as a standard or control in most immunoassay formats.

Disclaimer: For in vitro Laboratory Use Only. Not for diagnostic or therapeutic use. Not for human or animal consumption. Suggested applications of our products are not recommendations to use our products in violation of any patent or as a license under any patent of ImmunoReagents, Inc.

Product may not be resold or modified for resale without prior written approval of ImmunoReagents, Inc.